


Surgical Laser Workstations

More patients.
More treatments.
Less time.

Precisely.





The **next best thing**
in your practice.

When you are ready to expand treatment opportunities and your practice, Cynosure has the right surgical laser workstation for you. Large area to small area, laser lipolysis to cellulite reduction,¹ you can choose the platform that does what you need for the patients you want to attract.

Which Cynosure surgical laser is right for your practice?

Minimally Invasive Surgical Laser Treatment	Cellulaze	Smartlipo TriPlex	Smartlipo TriPlex + Cellulaze	PrecisionTx	Cellulaze + PrecisionTx	Smartlipo TriPlex + PrecisionTx	Smartlipo TriPlex + Cellulaze + PrecisionTx
Soft tissue contouring, such as neck				●	●	●	●
Soft tissue contouring, such as jawline and jowls				●	●	●	●
Subdermal ablation of glands, such as axillary sweat glands				●	●	●	●
Cellulite reduction ²	●		●		●		●
Laser-assisted lipolysis		●	●			●	●
Skin tightening ²		●	●			●	●
High-definition treatments		●	●			●	●

*It's been over a year and still
looks great...and feels great!*

– Debra R., Patient

Cellulaze™ Cellulite Laser Workstation

Transforming the way cellulite is treated

Cellulaze introduces a dramatic new way to treat cellulite by modifying its very structure with a single treatment. With its 1440nm Nd:YAG wavelength and a proprietary SideLaze3D™ side-firing laser fiber, ThermaGuide™ and SmartSense™ Cellulaze allows you to:⁸

- Thermally subcise the fibrous septae
- Target pockets of bulging fat
- Thicken and elasticize the dermis

Clinically proven

Cellulaze has been thoroughly evaluated to ensure safer and more effective results. Clinical study data⁸ showed:

- A single treatment improved the appearance of cellulite for one year with minimal side effects
- Increase in skin thickness by 25%
- Increase in skin elasticity by 29%
- Improvement in skin texture
- 93% patient satisfaction, and all would recommend Cellulaze treatment

Wavelength (Nd:YAG)
1440nm



Cellulaze gives women hope.

– Christine Petti, MD, Torrance, CA

Smartlipo TriPlex™

LaserBodySculpting™ Workstation

The gold standard in laser lipolysis

Smartlipo TriPlex has redefined laser lipolysis. Its three versatile wavelengths—1064nm, 1320nm and 1440nm—can be used independently or can be uniquely blended using Cynosure's MultiPlex™ technology to give you more treatment options.

The 1440nm Nd:YAG wavelength has 40 times greater absorption in fatty tissue than 924nm and 980nm diode lasers on the market.⁹

Smartlipo TriPlex is the first multiplex triple wavelength laser lipolysis workstation that's designed to treat all your patients' body contouring concerns:

- More effective laser lipolysis
- Tissue tightening achieved through heat-stimulated collagen generation²
- Reliable, controlled energy delivery
- Reduced treatment time and less downtime compared with traditional liposuction
- Less bleeding and bruising compared with traditional liposuction
- High-definition sculpting capabilities

Wavelengths (Nd:YAG)
1440nm
1064nm
1320nm



Smartlipo TriPlex has changed everything. It is predictable, fast, safe and more reliable. It causes no thermal damage, much less swelling and very little bruising. This is the standard of care in laser lipolysis.

– John Millard, MD, Millard Plastic Surgery, Lone Tree, CO

I finally have something to offer my patients who are not ready for a face-lift. The results have been fantastic, and all of my patients are thrilled.

– Dr. Beverly Fischer, Advanced Center for Plastic Surgery, Timonium, MD

PrecisionTx™

Precision Laser Workstation

Giving you new treatment opportunities from head-to-toe.


PrecisionTx from Cynosure introduces new minimally invasive laser procedures to your practice, enabling smaller area contouring such as jowls, jawline and neck. PrecisionTx treatment is ideal for patients who are looking for the best possible results with the least amount of downtime.³

One-time, one-hour treatments

PrecisionTx combines the Nd:YAG 1440nm wavelength with Cynosure's proprietary SideLaze800™ side-firing fiber and ThermoGuide™ technology to deliver laser energy subcutaneously to smaller areas. One-time, one-hour, minimally invasive office-based procedures provide desired outcomes—both to your patients and to your practice.

- Smaller area, head-to-toe precision treatments beyond laser lipolysis or cellulite reduction¹
- More convenient for smaller areas such as neck and jawline
- Subdermal ablation of glands, such as axillary sweat glands
- Proven results⁴ with minimal downtime⁵

Expand your surgical laser practice by giving it the PrecisionTx treatment.



Only Cynosure gives you so many choices.

Cynosure is committed to giving you more minimally invasive treatment opportunities than any other laser company. Leading laser-assisted lipolysis with Smartlipo TriPlex. The world's first longer-lasting one-time treatment for the appearance of cellulite³ with Cellulaze. And now, more accuracy for smaller-area laser treatments with PrecisionTx.

Contact your Cynosure Surgical Laser Consultant.

Customize a laser for your specific needs to optimize the potential of your practice.

**Call 800.886.2966
or visit www.cynosure.com**

References:

1. Cellulaze is intended for short-term improvement in the appearance of cellulite.
2. Tissue tightening achieved through heat-stimulated collagen generation. Orringer JS, Rittié L, Baker D, Voorhees JJ, Fisher G. Molecular mechanisms of nonablative fractionated laser resurfacing. *Br J Dermatol*. 2010 Oct;163(4):757-768
3. Gentile, R. Laser-assisted neck lift: high tech contouring and tightening. *Facial Plast Surg*. 2011 Aug 27 (4):331-345.
4. Sasaki, GH. "Clinical experience with the SideLaze800 fiber and the 1440nm wavelength Smartlipo TriPlex workstation for facial contouring," 2012 white paper. Twelve-patient study. Data on file.
5. Caplin, D. "Clinical evaluation and quantitative analysis of axillary hyperhidrosis treated with a unique targeted laser energy delivery method with 1 year follow up," ASLMS 2013 poster. Fifteen-patient study. Data on file.
6. Estimated absorption based on main chromophores in tissue: Dermis contains approximately 70% water. Human fatty adipose tissue contains approximately 75% lipid and 20% water. (Duck, FA. *Physical Properties of Tissue*. London: Academic Press; 1990; pp. 320-328.)
7. Levi JR, Veerappan A, Chen B, Mirkor M, Sierra R, Spiegel JH. Histologic evaluation of laser lipolysis comparing continuous wave vs pulsed lasers. *Arch Facial Plast Surg*. 2011;13(1):41-50.
8. DiBernardo, MD, Barry E. Treatment of cellulite using a 1440-nm pulsed laser with one-year follow-up. *Aesthet Surg J*. 2011 Mar 31 (3):328-341.
9. Estimation based on main chromophores in tissue; human fatty adipose tissue contains approximately 75% lipid and 20% water. (Duck, FA. *Physical Properties of Tissue*. London: Academic Press; 1990; p. 320.)



www.cynosure.com

CORPORATE HEADQUARTERS

United States

Cynosure, Inc.
5 Carlisle Road
Westford, MA 01886 USA
Tel: 978.256.4200
Toll-Free: 800.886.2966

